MEDUMAT Standard²

Ventilator

Safety information

Supplement to the instructions for use

This is a supplement to the instructions for use. Please read this safety information along with the instructions for use before using the product. Failure to observe the information contained therein can result in serious injuries or even death.

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Safety information

This safety information supplements the “Safety information” chapter of the instructions for use for the device. Please also always observe the instructions for use. The instructions for use and this supplement can be found at: www.weinmann-emergency.com.

WARNING

Risk of compromised therapy or therapy failure due to leaks or disconnections during ventilation!

Leaks during ventilation can reduce the applied tidal volume or airway pressure. If a serious leak occurs or there is a disconnection between the patient and the device, the patient cannot be ventilated. This may result in serious or life-threatening injury to the patient.

- ⇒ Use suitable monitoring (etCO₂, SpO₂ or expiratory volume measurement MVe).
- ⇒ Respond immediately to alarms and malfunctions that indicate leaks or disconnections (Airway pressure ↓, MVe ↓, Leakage on patient side (in CCSV mode)).
- ⇒ Check that the mask or airway device is correctly positioned.
- ⇒ Check all breathing circuit connections between the device and the patient for disconnections.

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Medical Technology GmbH