

COMBIBAG

One Resuscitator for Both, Adult and Child



☞ Make a 180-Degree Turn with COMBIBAG to Save Time, Space and Money



Thanks to the unique technology of COMBIBAG, the 2-in-1 resuscitator for adults and children, you always have the right equipment on hand. With a simple 180-degree turn of the COMBIBAG, you vary the volumes quickly and automatically adapt the tidal volume to the patient.

The safety valve limits the ventilation pressure and protects the patient during ventilation with a mask or tube. With COMBIBAG you save time, free up space in your emergency case and save on costs.

Your Benefits at a Glance

- Greater safety for users and patients: safety valve with two pressure levels for mask and tube ventilation
- One simple turn allows you to ventilate children from 10 kg using small volumes or adults using larger volumes
- Variety of connection options: improved oxygenation by means of oxygen reservoir or demand valve connection
- Simple cleaning and disinfection: single layer construction for optimum hygienic conditioning
- Maintenance and disassembly require no tools
- Latex-free: all components made of silicone or polysulfone

Product Versions



COMBIBAG
Resuscitator without mask
WM 11090



COMBIBAG
Resuscitator with size 5 mask
WM 11020

Accessories and Replacement Parts



1. Ventilation mask with self-filling cushion of silicone
 - For adults, size 5 WM 5074
 - For children and juveniles, size 3 WM 5082
2. OXYMAND demand valve with adapter and pressure hose 1.5 m
 - With connection nozzle G 3/8" and adapter for COMBIBAG WM 22127
 - With Walther plug and adapter for COMBIBAG WM 22128
 - With DIN 13260 plug and adapter for COMBIBAG WM 22129
3. PEEP valve, infinitely adjustable from 0 to 10 mbar WM 3215
4. Oxygen reservoir WM 11052

Technical Data

Product class as per directive 93 / 42 / EEC:	IIa
Area of application	
• Child segment:	10–16 kg body weight
• Adult segment:	> 16 kg body weight
Volume	
• Child segment:	up to 500 ml
• Adult segment:	500 to 1200 ml
Pressure limit:	switchable, 20 mbar or 60 mbar

Gas inlet connection for O ₂ supply:	Ø 6 mm
Patient connection:	15 mm / 22 mm
Dimensions (L x D):	approx. 330 mm x 125 mm
Weight:	approx. 410 g
Dead space:	5 ml
Materials	
• Valve:	Polysulfone
• Squeeze bag and valve membrane:	Silicone



Simply Professional

WEINMANN Emergency is a family-owned, internationally active medical technology company. With our mobile system solutions for emergency, transport and disaster medicine, we set standards for saving human lives. In close collaboration with professional users in emergency medical services, hospitals and military medical corps, we develop innovative medical products for ventilation and defibrillation. For more than 100 years we have offered our customers a high degree of reliability, extensive experience and quality made in Germany.

Headquarters

WEINMANN Emergency
Medical Technology GmbH + Co. KG
Frohbösestraße 12
22525 Hamburg
Germany

T: +49 40 88 18 96-0
F: +49 40 88 18 96-480
T: +49 40 88 18 96-120 Customer Service
T: +49 40 88 18 96-122 After-Sales Service
E: info@weinmann-emt.de

Center for Production, Logistics, Service

WEINMANN Emergency
Medical Technology GmbH + Co. KG
Siebenstücken 14
24558 Henstedt-Ulzburg
Germany

China

Weinmann (Shanghai) Medical Device Trading Co. Ltd.
T: +86 21 52 30 22 25 • info@weinmann-emt.cn

U.A.E.

WEINMANN Emergency Medical Technology GmbH + Co.KG (Branch)
T: +971 432 100 31 • info-dubai@weinmann-emt.com

France

WEINMANN Emergency France SARL – Paris – Les Ulis
T: +33 1 69 41 51 20 • info@weinmann-emt.fr

Russia

Weinmann SPb GmbH – St. Petersburg
T: +7 812 633 30 82 • info@weinmann-emt.ru

Singapore

Weinmann Singapur PTE, Ltd.
T: +65 65 09 44 30 • info-singapore@weinmann-emt.sg

Spain

WEINMANN Emergency Medical Technology GmbH + Co. KG
T: +34 91 79 01 137 • info-spain@weinmann-emt.es

USA

Weinmann Emergency LP
T: +1 770-274-2417 • info@weinmann-emergency.com