



Price list

Valid from 01.04.2018

Non-binding, recommended sales prices. All prices excluding sales tax.

This price list remains valid until a subsequent price list is issued. Old price lists lose their validity on the date the latest price list goes into effect.

Customized Configurations

This price list may contain sample customized configurations, which are identified by the addition of a four-digit configuration number to the article number. In the example WM XXXXX-YYYY, the last four digits representing the customized configuration are separated from the actual article number by a hyphen.

All the prices for the configurations are based on a sample configuration of our choice generally numbered -0000). Your final price can differ, depending on your customized configuration!

Articles which are not configured will still be presented this way:

WM XXXXX (without a hyphen).

Up-to-date with our Customer Login

Our price lists are regularly updated at short notice. Before you place an order, please obtain information about the currently valid price lists by signing into the Customer Login at www.weinmann-emergency.de or by calling our Customer Service.

Shipping + Ordering of Minimum Quantities

A minimum quantity surcharge of 30 € applies to deliveries in foreign countries with a total value of less than 250 €.

All rights reserved

We reserve the right to make changes to constructions, designs, equipment and prices, particularly during ongoing development of new production series. Without our express permission, any form of exploitation and use of copyright for and from this publication is restricted solely to WEINMANN Emergency Medical Technology GmbH + Co. KG.

Do you have any questions?

Our Customer Service staff will be happy to help you.

Germany

WEINMANN Emergency
Medical Technology GmbH + Co. KG
Frohbösestraße 12
22525 Hamburg
Germany

T: +49 40 88 18 96-0
F: +49 40 88 18 96-480
T: +49 40 88 18 96-120 Customer Service
T: +49 40 88 18 96-122 After-Sales Service
E: info@weinmann-emt.de

Center for Production, Logistics, Service
WEINMANN Emergency
Medical Technology GmbH + Co. KG
Siebenstücken 14
24558 Henstedt-Ulzburg
Germany

China

Weinmann (Shanghai)
Medical Device Trading Co. Ltd.
T: +86 21 52 30 22 25
E: info@weinmann-emt.cn

U.A.E.

WEINMANN Emergency
Medical Technology GmbH + Co.KG (Branch)
T: +971 (0) 432 100 31
E: info-dubai@weinmann-emt.com

Frankreich

WEINMANN Emergency France SARL
T: +33 1 69 41 51 20
E: info@weinmann-emt.fr

Russia

Weinmann SPb GmbH – St. Petersburg
T: +7 812 633 30 82
E: info@weinmann-emt.ru

Singapore

Weinmann Singapur PTE, Ltd.
T: +65 65 09 44 30
E: info-singapore@weinmann-emt.sg

Spain

WEINMANN Emergency
Medical Technology GmbH + Co. KG
T: +34 91 79 01 137
E: info-spain@weinmann-emt.es



Table of content

Oxygen Systems for Emergency Vehicles	1
Oxygen Pressure Reducer Systems	1
OXYWAY Fix I	1
OXYWAY Fix III	2
OXYWAY Fine I	6
OXYWAY Fast I	9
OXYWAY Fast II	11
OXYWAY Fast III	16
OXYWAY Flow Dosage	17
OXYWAY Click	17
Accessories and Replacement Parts for all pressure reducers	17
Consumables for inhalation	18
Oxygen Control Panels	18
Oxygen administrators	18
Accessories and Replacement Parts for Oxygen Control Panels	19
O2 plugs and adapters	19
Connection Nozzles	21
Distributors	21
Pressure hoses	22
Oxygen Cylinders	25
Oxygen Cylinders 0.3 L	25
Oxygen Cylinders 0.8 L	26
Oxygen Cylinders 2 L	26
Oxygen Cylinders 2 L Aluminum	27
Oxygen Cylinders 5 L	27
Oxygen Cylinders 10 L	27
Accessories for Oxygen Cylinders	28
O2 Box Rescue	29
Accessories and Replacement Parts	31

Oxygen Systems for Emergency Vehicles

Oxygen Pressure Reducer Systems

*Further connection threads
are available upon request*

OXYWAY Fix I

Connection bolt thread: G 3/4

- WM 30101-0000 OXYWAY Fix I, G 3/4
Oxygen pressure reducer
- High-pressure manual connection G 3/4
 - 30 mm connection bolt, right-hand side
 - pressure gauge
 - permanently set to 4.5 bar
 - Flow outlet (4 l/min), bottom, UNF 9/16

Connection bolt thread: G 5/8, 60° Bullnose

- WM 30102-0000 OXYWAY Fix I, G 5/8, 60° Bullnose
Oxygen pressure reducer
- High-pressure manual connection G 5/8, 60° Bullnose
 - 30 mm connection bolt, right-hand side
 - with pressure gauge
 - permanently set to 4.5 bar
 - Flow outlet (4 l/min), bottom, UNF 9/16

Connection bolt thread: G 5/8, 72° Bullnose

- WM 30103-0000 OXYWAY Fix I, G 5/8, 72° Bullnose
Oxygen pressure reducer
- High-pressure manual connection G 5/8, 72° Bullnose
 - 30 mm connection bolt, right-hand side
 - with pressure gauge
 - permanently set to 4.5 bar
 - Flow outlet (4 l/min), bottom, UNF 9/16

- WM 30123-0000 OXYWAY Fix I, G 5/8, 72° Bullnose
Oxygen pressure reducer
- High-pressure manual connection G 5/8, 72° Bullnose
 - 50 mm connection bolt, right-hand side
 - with pressure gauge
 - permanently set to 4.5 bar

Connection bolt thread: W21,8 x 1/14"

- WM 30104-0000 "OXYWAY Fix I, W21,8 x 1/14"
Oxygen pressure reducer
- High-pressure manual connection W21,8 x 1/14"
 - 30 mm connection bolt, right-hand side
 - with pressure gauge
 - permanently set to 4.5 bar
 - Flow outlet (4 l/min), bottom, UNF 9/16

Connection bolt thread: NGO

- WM 30115-0000 OXYWAY Fix I, CGA 540, NGO
Oxygen pressure reducer
- High-pressure manual connection CGA 540, NGO
 - 40 mm connection bolt, right-hand side
 - with pressure gauge
 - permanently set to 4.5 bar
 - Flow outlet (4 l/min), bottom

Connection bolt thread: CGA 870 PIN-Index

- WM 30106-0000 OXYWAY Fix I, PIN-Index
Oxygen pressure reducer
- High-pressure manual PIN-Index connection
 - 30 mm connection bolt, right-hand side
 - with pressure gauge
 - permanently set to 4.5 bar
 - Flow outlet (4 l/min), bottom, UNF 9/16

Connection bolt thread: G 5/8, 90° Bullnose

- WM 30107-0000 OXYWAY Fix I, G 5/8, 90° Bullnose
Oxygen pressure reducer
- High-pressure manual connection G 5/8, 90° Bullnose (BE/LU)
 - 30 mm connection bolt, right-hand side
 - with pressure gauge
 - permanently set to 4.5 bar

OXYWAY Fix III

for MEDUMAT / MODUL-Systems:

Connection bolt thread: G 3/4

- WM 30301-0000 OXYWAY Fix III, G 3/4
Oxygen pressure reducer
- High-pressure G 3/4 manual connection
 - 30 mm connection bolt, right-hand side
 - Pressure gauge
 - set to 4.5 bar,
 - G 3/8 pressure outlet (190 l/min), bottom
- WM 30381-0000 OXYWAY Fix III, G 3/4 with flat contents gauge
Pressure reducer for oxygen,
- High-pressure manual connection G 3/4
 - 20 mm short connection bolts, right
 - Contents gauge (flat, zero flow position at 12 o'clock)
 - Permanently set at 4.5 bar
 - Flow outlet G 3/8 (190 l/min), below
- WM 30391-0000 OXYWAY Fix III, G 3/4 with flat contents gauge
Pressure reducer for oxygen,
- High-pressure manual connection G 3/4
 - 20 mm short connection bolts, left
 - Contents gauge (flat, zero flow position at 6 o'clock)
 - Permanently set at 4.5 bar
 - Flow outlet G 3/8 (190 l/min), below
- WM 30311-0000 OXYWAY Fix III, G 3/4 with flat gauge
Oxygen pressure reducer,
- High-pressure G 3/4 manual connection
 - 30 mm connection bolt, right-hand side
 - Pressure gauge (flat, zero-position 3 o'clock)
 - set to 4.5 bar,
 - G 3/8 pressure outlet (190 l/min), bottom
- WM 30361-0000 OXYWAY Fix III, G 3/4
Oxygen pressure reducer
- High-pressure G3/4 manual connection
 - 50 mm connection bolt, left-hand side
 - Pressure gauge
 - set to 4.5 bar,
 - G 3/8 pressure outlet (190 l/min), bottom
- WM 30341-0000 OXYWAY Fix III, G 3/4 with flat gauge
Oxygen pressure reducer,
- High-pressure G 3/4 manual connection
 - 30 mm connection bolt, left-hand side

	<ul style="list-style-type: none"> - Pressure gauge (flat, zero-position 12 o'clock) - set to 4.5 bar, - G 3/8 pressure outlet (190 l/min), bottom
WM 30351-0000	<p>OXYWAY Fix III, G 3/4 Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure G 3/4 manual connection - 30 mm connection bolt, left-hand side - Pressure gauge - set to 4.5 bar, - G 3/8 pressure outlet (190 l/min), bottom
WM 30321-0000	<p>OXYWAY Fix III, G 3/4 Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure G 3/4 manual connection - 50 mm connection bolt, right-hand side - Pressure gauge - set to 4.5 bar, - G 3/8 pressure outlet (190 l/min), bottom
	<p>Connection bolt thread: G 5/8, 60° Bullnose</p>
WM 30302-0000	<p>OXYWAY Fix III, G 5/8, 60° Bullnose Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G5/8, 60° Bullnose - 30 mm connection bolt, right-hand side, - Pressure gauge - set to 4.5 bar,
WM 30342-0000	<p>OXYWAY Fix III, G 5/8, 60° Bullnose, with flat gauge Oxygen pressure reducer,</p> <ul style="list-style-type: none"> - High-pressure manual connection G 5/8, 60° Bullnose - 30 mm connection bolt, left-hand side - Pressure gauge (flat, zero-position 12 o'clock) - set to 4.5 bar, - G 3/8 pressure outlet (190 l/min), bottom
WM 30322-0000	<p>OXYWAY Fix III, G 5/8, 60° Bullnose Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G 5/8, 60° Bullnose - 50 mm connection bolt, right-hand side, - Pressure gauge - set to 4.5 bar, - G 3/8 pressure outlet (190 l/min), bottom
	<p>Connection bolt thread: G 5/8, 72° Bullnose</p>
WM 30303-0000	<p>OXYWAY Fix III, G 5/8, 72° Bullnose Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G 5/8, 72° Bullnose - 30 mm connection bolt, right-hand side, - Pressure gauge - set to 4.5 bar, - G 3/8 pressure outlet (190 l/min), bottom
WM 30343-0000	<p>OXYWAY Fix III, G 5/8, 72° Bullnose, with flat gauge Oxygen pressure reducer,</p> <ul style="list-style-type: none"> - High-pressure manual connection G 5/8, 72° Bullnose - 30 mm connection bolt, left-hand side - Pressure gauge (flat, zero-position 12 o'clock) - set to 4.5 bar,

	- G 3/8 pressure outlet (190 l/min), bottom
WM 30323-0000	OXYWAY Fix III, G 5/8, 72° Bullnose Oxygen pressure reducer - High-pressure manual connection G 5/8, 72° Bullnose - 50 mm connection bolt, right-hand side, - Pressure gauge - set to 4.5 bar, - G 3/8 pressure outlet (190 l/min), bottom
WM 30363-0000	OXYWAY Fix III, G 5/8, 72° Bullnose Oxygen pressure reducer - High-pressure manual connection G 5/8, 72° Bullnose - 50 mm connection bolt, left-hand side, - Pressure gauge - set to 4.5 bar, - G 3/8 pressure outlet (190 l/min), bottom
	Connection bolt thread: W21,8 x 1/14"
WM 30304-0000	"OXYWAY Fix III, W21,8 x 1/14"" Oxygen pressure reducer - High-pressure manual connection W21,8 x 1/14"" - 30 mm connection bolt, right-hand side, - Pressure gauge - set to 4.5 bar, - G 3/8 pressure outlet (190 l/min), bottom "
WM 30354-0000	"OXYWAY Fix III, W21,8 x 1/14"" Oxygen pressure reducer - High-pressure manual connection W21,8 x 1/14"" - 30 mm connection bolt, left-hand side, - Pressure gauge - set to 4.5 bar, - G 3/8 pressure outlet (190 l/min), bottom "
WM 30344-0000	"OXYWAY Fix III, W21,8 x 1/14""", with flat gauge Oxygen pressure reducer, - High-pressure manual connection W21,8 x 1/14"" - 30 mm connection bolt, left-hand side - Pressure gauge (flat, zero-position 12 o'clock) - set to 4.5 bar, - G 3/8 pressure outlet (190 l/min), bottom "
WM 30324-0000	"OXYWAY Fix III, W21,8 x 1/14"" Oxygen pressure reducer - High-pressure manual connection W21,8 x 1/14"" - 50 mm connection bolt, right-hand side, - Pressure gauge - set to 4.5 bar, - G 3/8 pressure outlet (190 l/min), bottom "
WM 30364-0000	"OXYWAY Fix III, W21,8 x 1/14"" Oxygen pressure reducer - High-pressure manual connection W21,8 x 1/14"" - 50 mm connection bolt, left-hand side, - Pressure gauge - set to 4.5 bar, - G 3/8 pressure outlet (190 l/min), bottom "
	Connection bolt thread: NGO
WM 30305-0000	OXYWAY Fix III, CGA 540, NGO

Oxygen pressure reducer
 - High-pressure manual connection CGA 540, NGO
 - 30 mm connection bolt, right-hand side,
 - Pressure gauge
 - set to 4.5 bar,
 - G 3/8 pressure outlet (190 l/min), bottom

WM 30375-0000 OXYWAY Fix III, CGA 540, NGO, with flat gauge
 Oxygen pressure reducer,
 - High-pressure manual connection CGA 540, NGO
 - 40 mm connection bolt, left-hand side
 - Pressure gauge (flat, zero-position 12 o'clock)
 - set to 4.5 bar,
 - G 3/8 pressure outlet (190 l/min), bottom

WM 30325-0000 OXYWAY Fix III, CGA 540, NGO
 Oxygen pressure reducer
 - High-pressure manual connection CGA 540, NGO
 - 50 mm connection bolt, right-hand side,
 - Pressure gauge
 - set to 4.5 bar,
 - G 3/8 pressure outlet (190 l/min), bottom

Connection bolt thread: CGA 870 PIN-Index

WM 30306-0000 OXYWAY Fix III, PIN-Index
 Oxygen pressure reducer
 - High-pressure manual PIN-Index connection
 - 30 mm connection bolt, right-hand side,
 - Pressure gauge
 - set to 4.5 bar,
 - G 3/8 pressure outlet (190 l/min), bottom

WM 30346-0000 OXYWAY Fix III, PIN-Index, with flat gauge
 Oxygen pressure reducer,
 - High-pressure PIN-Index manual connection
 - 30 mm connection bolt, left-hand side
 - Pressure gauge (flat, zero-position 12 o'clock)
 - set to 4.5 bar,
 - G 3/8 pressure outlet (190 l/min), bottom

WM 30326-0000 OXYWAY Fix III, PIN-Index
 Oxygen pressure reducer
 - High-pressure manual PIN-Index connection
 - 50 mm connection bolt, right-hand side,
 - Pressure gauge
 - set to 4.5 bar,
 - G 3/8 pressure outlet (190 l/min), bottom

Connection bolt thread: G 5/8, 90° Bullnose

WM 30307-0000 OXYWAY Fix III, G 5/8, 90° Bullnose
 Oxygen pressure reducer
 - High-pressure manual connection, G 5/8,
 90° Bullnose
 - 30 mm connection bolt, right-hand side
 - Pressure gauge
 - set to 4.5 bar,
 - G 3/8 pressure outlet (190 l/min), bottom

WM 30347-0000 OXYWAY Fix III, G 5/8, 90° Bullnose, with flat
 gauge
 Oxygen pressure reducer,
 - High-pressure manual connection G 5/8,
 90° Bullnose
 - 30 mm connection bolt, left-hand side
 - Pressure gauge (flat, zero-position 12 o'clock)
 - set to 4.5 bar,
 - G 3/8 pressure outlet (190 l/min), bottom

Connection bolt thread: NF E29-650/F 90° Bullnose

WM 30308-0000 OXYWAY Fix III, NF E29-650/F (AFNOR)
Oxygen pressure reducer
- High-pressure manual connection NF E29-650/F (AFNOR)
- 30 mm connection bolt, right-hand side
- Pressure gauge
- set to 4,5 bar,
- G 3/8 pressure outlet (190 l/min), bottom

WM 30421-0000 OXYWAY Fix III, G 3/4
Oxygen pressure reducer
- High-pressure manual connection G 3/4
- 50 mm connection bolt, right-hand side
- Pressure gauge
- permanently set to 4,5 bar
- Pressure outlet: DIN coupling (150 l/min), bottom

Connection bolt thread: G 5/8, 72° Bullnose

WM 30413-0000 OXYWAY Fix III, G 5/8, 72° Bullnose
Oxygen pressure reducer
- High-pressure manual connection G 5/8, 72° Bullnose
- 30 mm connection bolt, right-hand side
- Pressure gauge
- permanently set to 4,5 bar
- Pressure outlet: bayonet coupling (120 l/min), bottom

WM 30423-0000 OXYWAY Fix III, G 5/8, 72° Bullnose
Oxygen pressure reducer
- High-pressure manual connection G 5/8, 72° Bullnose
- 50 mm connection bolt, right-hand side
- Pressure gauge
- permanently set to 4,5 bar
- Pressure outlet: bayonet coupling (120 l/min), bottom

Connection bolt thread: CGA 870 PIN-Index

WM 30416-0000 OXYWAY Fix III, PIN-Index
Oxygen pressure reducer
- High-pressure manual PIN-Index connection
- 30 mm connection bolt, right-hand side
- Pressure gauge
- permanently set to 4,5 bar
- Pressure outlet: bayonet coupling (120 l/min), bottom

Connection bolt thread: G 5/8, 90° Bullnose

WM 30417-0000 OXYWAY Fix III, G 5/8, 90° Bullnose
Oxygen pressure reducer
- High-pressure manual connection G 5/8, 90° Bullnose
- 30 mm connection bolt, right-hand side
- Pressure gauge
- permanently set to 4,5 bar
- Pressure outlet: DIN coupling (150 l/min), bottom

OXYWAY Fine I**Connection bolt thread: G 3/4**

WM 30501-0000 OXYWAY Fine I, G 3/4
Oxygen pressure reducer

	<ul style="list-style-type: none"> - High-pressure manual connection G 3/4 - 30 mm connection bolt, right-hand side - Contents and flow gauge - Infinitely adjustable - Operating range 3 - 15 L/min (0.5 - 4.5 bar) - flow outlet (3-15l/min), bottom, UNF 9/16
WM 30521-0000	<p>OXYWAY Fine I, G 3/4 Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G 3/4 - 50 mm connection bolt, right-hand side - Contents and flow gauge - Infinitely adjustable - Operating range 3 - 15 L/min (0.5 - 4.5 bar) - flow outlet (3-15l/min), bottom, UNF 9/16
WM 30531-0000	<p>OXYWAY Fine I, G 3/4 Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G 3/4 - 80 mm connection bolt, right-hand side - Contents and flow gauge - Infinitely adjustable - Operating range 3 - 15 L/min (0.5 - 4.5 bar) - flow outlet (3-15l/min), bottom, UNF 9/16
WM 30541-0000	<p>OXYWAY Fine I Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G 3/4 - 100 mm connection bolt, right-hand side - Contents and flow gauge - Infinitely adjustable - Operating range 3 - 15 L/min (0.5 - 4.5 bar) - flow outlet (3-15l/min), bottom, outlet UNF 9/16
	<p>Connection bolt thread: G 5/8, 60° Bullnose</p>
WM 30502-0000	<p>OXYWAY Fine I, G 5/8, 60° Bullnose pressure reducer for oxygen</p> <ul style="list-style-type: none"> - High pressure manual connection G 5/8, 60° Bullnose - 30 mm connection bolt, right - Content and flow gauge - infinitely adjustable - Operating range 3 - 15 l/min (0.5 - 4.5 bar) - Flow outlet (3-15l/min), UNF 9/16
WM 30522-0000	<p>OXYWAY Fine I, G 5/8, 60° Bullnose Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G 5/8, 60° Bullnose - 50 mm connection bolt, right-hand side - Contents and flow gauge - Infinitely adjustable - Operating range 3 - 15 L/min (0.5 - 4,5 bar) - flow outlet (3-15l/min), bottom, UNF 9/16
	<p>Connection bolt thread: G 5/8, 72° Bullnose</p>
WM 30503-0000	<p>OXYWAY Fine I, G 5/8, 72° Bullnose Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G 5/8, 72° Bullnose - 30 mm connection bolt, right-hand side - Contents and flow gauge - Infinitely adjustable - Operating range 3 - 15 L/min (0.5 - 4,5 bar) - flow outlet (3-15l/min), bottom, UNF 9/16
WM 30523-0000	<p>OXYWAY Fine I, G 5/8, 72° Bullnose</p>

- Oxygen pressure reducer
- High-pressure manual connection G 5/8, 72° Bullnose
 - 50 mm connection bolt, right-hand side
 - Contents and flow gauge
 - Infinitely adjustable
 - Operating range 3 - 15 L/min (0.5 - 4,5 bar)
 - flow outlet (3-15l/min), bottom, UNF 9/16

Connection bolt thread: W21,8 x 1/14"

- WM 30504-0000 "OXYWAY Fine I, W21,8 x 1/14""
- Oxygen pressure reducer
- High-pressure manual connection W21,8 x 1/14""
 - 30 mm connection bolt, right-hand side
 - Contents and flow gauge
 - Infinitely adjustable
 - Operating range 3 - 15 L/min (0.5 - 4,5 bar)
 - flow outlet (3-15l/min), bottom, UNF 9/16

- WM 30524-0000 "OXYWAY Fine I, W21,8 x 1/14""
- Oxygen pressure reducer
- High-pressure manual connection W21,8 x 1/14""
 - 50 mm connection bolt, right-hand side
 - Contents and flow gauge
 - Infinitely adjustable
 - Operating range 3 - 15 L/min (0.5 - 4,5 bar)
 - flow outlet (3-15l/min), bottom, UNF 9/16

Connection bolt thread: NGO

- WM 30515-0000 OXYWAY Fine I, CGA 540, NGO
- Oxygen pressure reducer
- High-pressure manual connection CGA 540, NGO
 - 30 mm connection bolt, right-hand side
 - Contents and flow gauge
 - Infinitely adjustable
 - Operating range 3 - 15 L/min (0.5 - 4,5 bar)
 - flow outlet (3-15l/min), bottom, UNF 9/16

Connection bolt thread: CGA 870 PIN-Index

- WM 30506-0000 OXYWAY Fine I, PIN-Index
- Oxygen pressure reducer
- High-pressure Pin Index connection
 - 30 mm connection bolt, right-hand side
 - Contents and flow gauge
 - Infinitely adjustable
 - Operating range 3 - 15 L/min (0.5 - 4,5 bar)
 - flow outlet (3-15l/min), bottom, UNF 9/16

- WM 30526-0000 OXYWAY Fine I, PIN-Index
- Oxygen pressure reducer
- High-pressure Pin Index connection
 - 50 mm connection bolt, right-hand side
 - Contents and flow gauge
 - Infinitely adjustable
 - Operating range 3 - 15 L/min (0.5 - 4,5 bar)
 - flow outlet (3-15l/min), bottom, UNF 9/16

Connection bolt thread: G 5/8, 90° Bullnose

- WM 30507-0000 OXYWAY Fine I, G 5/8, 90° Bullnose
- Oxygen pressure reducer
- High-pressure manual connection G 5/8, 90° Bullnose
 - 30 mm connection bolt, right-hand side

- Contents and flow gauge
- Infinitely adjustable
- Operating range 3 - 15 L/min (0.5 - 4,5 bar)
- flow outlet (3-15l/min), bottom, UNF 9/16

OXYWAY Fast I

Connection bolt thread: G 3/4

WM 30861-0000

- OXYWAY Fast I, G 3/4
Pressure reducer for oxygen
- High pressure manual connection G 3/4
 - 20 mm connection bolt, right
 - Dosage dial and contents gauge
 - Flow outlet (1-15l/min), bottom, UNF 9/16
 - dosages in levels 1,2,3,4,5,6,9,12,15 l/min

WM 30801-0000

- OXYWAY Fast I, G 3/4
Oxygen pressure reducer
- High-pressure manual connection G 3/4
 - 30 mm connection bolt, right-hand side
 - Dosage dial and pressure gauge
 - flow outlet (1-15l/min), bottom, UNF 9/16
 - Dosages in levels 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min

WM 30811-0000

- OXYWAY Fast I, G 3/4 with flat gauge,
Oxygen pressure reducer,
- High-pressure manual connection G 3/4
 - 50 mm connection bolt, right-hand side
 - Dosage dial and pressure gauge (flat, zero-position 3 o'clock)
 - flow outlet (1-15l/min), sideways, UNF 9/16
 - Dosages in levels 1,2,3,4,5,6, 9, 12, 15 l/min

WM 33801-0000

- OXYWAY Fast I, G 3/4 with flat gauge,
Oxygen pressure reducer,
- High-pressure manual connection G 3/4
 - 50 mm connection bolt, left-hand side
 - Dosage dial and pressure gauge (flat, zero-position 3 o'clock)
 - flow outlet (1-15l/min), up, UNF 9/16
 - Dosages in levels 1,2,3,4,5,6, 9, 12, 15 l/min

WM 33781-0000

- OXYWAY Fast I, G 3/4
Oxygen pressure reducer,
- High-pressure manual connection G 3/4
 - 30 mm connection bolt, left-hand side
 - Dosage dial and pressure gauge (flat, zero-position 9 o'clock)
 - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16, Outlet right, 270° turned

WM 33991-0000

- OXYWAY Fast I, G 3/4
Oxygen pressure reducer,
- High-pressure manual connection G 3/4
 - 50 mm connection bolt, left-hand side
 - Dosage dial and pressure gauge (flat, zero-position 9 o'clock)
 - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16

Connection bolt thread: G 5/8,60° Bullnose

WM 30802-0000

- OXYWAY Fast I, G 5/8, 60° Bullnose
Oxygen pressure reducer
- High pressure manual connection G 5/8, 60° Bullnose

- 30 mm connection bolt, right-hand side
- Dosage dial and pressure gauge
- flow outlet (1-15l/min), bottom, UNF 9/16
- Dosages in levels 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min

Connection bolt thread: G 5/8, 72° Bullnose

- WM 30803-0000 OXYWAY Fast I, G 5/8, 72° Bullnose
Oxygen pressure reducer
- High-pressure manual connection G 5/8, 72° Bullnose
 - 30 mm connection bolt, right-hand side
 - Dosage dial and pressure gauge
 - flow outlet (1-15l/min), bottom, UNF 9/16
 - Dosages in levels 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min

- WM 30823-0000 OXYWAY Fast I, G 5/8, 72° Bullnose
Oxygen pressure reducer
- High-pressure manual connection G 5/8, 72° Bullnose
 - 50 mm connection bolt, right-hand side
 - Dosage dial and pressure gauge
 - flow outlet (1-15l/min), bottom, UNF 9/16
 - Dosages in levels 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min

Connection bolt thread: W21,8 x 1/14"

- WM 30804-0000 "OXYWAY Fast I, W21,8 x 1/14"
Oxygen pressure reducer
- High-pressure manual connection W21,8 x 1/14"
 - 30 mm connection bolt, right-hand side
 - Dosage dial and pressure gauge
 - Flow outlet (1-15l/min), bottom, UNF 9/16
 - Dosages in levels 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min

- WM 30805-0000 OXYWAY Fast I, CGA 540, NGO
Oxygen Pressure Reducer
- High-pressure manual connection W21,8/21,7
 - 30 mm connection bolt, right-hand side
 - Dosage dial and pressure gauge
 - Flow outlet (1-15l/min), bottom, UNF 9/16
 - Dosages in levels 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min

- WM 30814-0000 OXYWAY Fast I, W21,8 x 1/14"
Oxygen pressure reducer
- High-pressure manual connection W21,8 x 1/14"
 - 30 mm connection bolt, right-hand side
 - Dosage dial and pressure gauge, flat
 - Flow outlet (1-15l/min), bottom, UNF 9/16
 - Dosages in levels 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min

Connection bolt thread: NGO

- WM 30815-0000 OXYWAY Fast I, CGA 540, NGO
Oxygen pressure reducer
- High-pressure manual connection CGA 540, NGO
 - 40 mm connection bolt, right-hand side
 - Dosage dial and pressure gauge
 - flow outlet (1-15l/min), bottom, UNF 9/16
 - Dosages in levels 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min

Connection bolt thread: CGA 870 PIN-Index

- WM 30806-0000 OXYWAY Fast I, PIN-Index
Oxygen pressure reducer

- High-pressure Pin Index connection
- 30 mm connection bolt, right-hand side
- Dosage dial and pressure gauge
- flow outlet (1-15l/min), bottom, UNF 9/16
- Dosages in levels 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min

Connection bolt thread: G 5/8, 90° Bullnose

- WM 30807-0000 OXYWAY Fast I, G 5/8, 90° Bullnose
Oxygen pressure reducer
- High-pressure manual connection G 5/8, 90° Bullnose
 - 30 mm connection bolt, right-hand side
 - Dosage dial and pressure gauge
 - flow outlet (1-15l/min), bottom, UNF 9/16
 - Dosages in levels 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min

OXYWAY Fast II

Connection bolt thread: G 3/4

- WM 33811-0000 OXYWAY Fast II, G 3/4
Oxygen Pressure Reducer
- High-pressure manual connection G 3/4
 - 20 mm connection bolt, right-hand side,
 - Adjustable dosage dial and pressure gauge
 - G 3/8 pressure outlet (190 l/min), bottom, with locking screw
 - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 l/min, UNF 9/16 turned 270°

- WM 33821-0000 OXYWAY Fast II, G 3/4
Oxygen pressure reducer
- High-pressure manual connection G 3/4
 - 20 mm connection bolt, left-hand side
 - Dosage dial and pressure gauge
 - G 3/8 pressure outlet(190 L/min), bottom
 - Flow outlet, dadjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16, turned 90°

- WM 31831-0000 OXYWAY Fast II, G 3/4
Pressure reducer for oxygen,
- High-pressure manual connection G 3/4
 - 20 mm short connection bolt, left
 - Dosage dial and contents gauge
 - Pressure outlet G 3/8, bottom, free flow (about 190 l/min), with locking screw
 - Flow outlet, dosages set at Levels 1, 2, 3, 4, 5, 6, 9, 12, 15 l/min, UNF 9/16, bottom

- WM 30851-0000 OXYWAY Fast II, G 3/4
Oxygen pressure reducer
- High-pressure manual connection G 3/4
 - 30 mm connection bolt, right-hand side,
 - Adjustable dosage dial and pressure gauge
 - G 3/8 pressure outlet (190 l/min), bottom, with locking screw
 - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 l/min, UNF 9/16

- WM 30871-0000 OXYWAY Fast II, G 3/4
Oxygen pressure reducer
- High-pressure manual connection G 3/4
 - 50 mm connection bolt, right-hand side,
 - Adjustable dosage dial and pressure gauge
 - G 3/8 pressure outlet (190 l/min), bottom, with locking screw
 - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 l/min, UNF 9/16

WM 30841-0000	<p>OXYWAY Fast II, G 3/4 Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G 3/4 - 30 mm connection bolt, right-hand side - Dosage dial and pressure gauge - Pressure outlet: DIN coupling (150 l/min), bottom - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 l/min, UNF 9/16
WM 30881-0000	<p>OXYWAY Fast II, G 3/4 Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G 3/4 - 50 mm connection bolt, right-hand side - Dosage dial and pressure gauge (flat, zero-position 3 o'clock) - Pressure outlet: DIN coupling (150 l/min), bottom - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 l/min, UNF 9/16, sideways, 90° turned
WM 30891-0000	<p>OXYWAY Fast II, G 3/4 Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G 3/4 - 30 mm connection bolt, right-hand side - Dosage dial and pressure gauge - G 3/8 pressure outlet(190 L/min), bottom - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16, turned 270°, right <p>Connection bolt thread: G 5/8,60° Bullnose</p>
WM 30852-0000	<p>OXYWAY Fast II, G 5/8, 60° Bullnose Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G 5/8, 60° Bullnose - 30 mm connection bolt, right-hand side, - Adjustable dosage dial and pressure gauge - G 3/8 pressure outlet (190 l/min), bottom, with locking screw - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 l/min, UNF 9/16
WM 30892-0000	<p>OXYWAY Fast II, G 5/8, 60° Bullnose Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G 5/8, 60° Bullnose - 30 mm connection bolt, right-hand side - Dosage dial and pressure gauge - G 3/8 pressure outlet (190 L/min), bottom - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16, turned 270°, right <p>Connection bolt thread: G 5/8, 72° Bullnose</p>
WM 30853-0000	<p>OXYWAY Fast II, G 5/8, 72° Bullnose Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G 5/8, 72° Bullnose - 30 mm connection bolt, right-hand side, - Adjustable dosage dial and pressure gauge - G 3/8 pressure outlet (190 l/min), bottom, with locking screw - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 l/min, UNF 9/16
WM 30873-0000	<p>OXYWAY Fast II, G 5/8, 72° Bullnose Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G 5/8, 72° Bullnose - 50 mm connection bolt, right-hand side,

	<ul style="list-style-type: none"> - Adjustable dosage dial and pressure gauge - G 3/8 pressure outlet (190 l/min), bottom, with locking screw - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 l/min, UNF 9/16
WM 30893-0000	<p>OXYWAY Fast II, G 5/8, 72° Bullnose Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G 5/8, 72° Bullnose - 30 mm connection bolt, right-hand side - Dosage dial and pressure gauge - G 3/8 pressure outlet (190 l/min), bottom - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16, turned 270°, right <p>Connection bolt thread: W21,8 x 1/14"</p>
WM 30854-0000	<p>"OXYWAY Fast II, W21,8 x 1/14"" Oxygen Pressure Reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection W21,8 x 1/14"" - 30 mm connection bolt, right-hand side, - Adjustable dosage dial and pressure gauge - G 3/8 pressure outlet (190 l/min), bottom, with locking screw - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 l/min, UNF 9/16
WM 30894-0000	<p>"OXYWAY Fast II, W21,8 x 1/14"" Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection W21,8 x 1/14"" - 30 mm connection bolt, right-hand side - Dosage dial and pressure gauge - G 3/8 pressure outlet (190 l/min), bottom - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16" turned 270°, right <p>Connection bolt thread: NGO</p>
WM 30865-0000	<p>OXYWAY Fast II, CGA 540, NGO Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection CGA 540, NGO - 30 mm connection bolt, right-hand side, - Adjustable dosage dial and pressure gauge - G 3/8 pressure outlet (190 l/min), bottom, with locking screw - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 l/min, UNF 9/16 <p>Connection bolt thread: CGA 870 PIN-Index</p>
WM 30856-0000	<p>OXYWAY Fast II, PIN-Index Oxygen Pressure Reducer</p> <ul style="list-style-type: none"> - High-pressure PIN-Index manual connection - 30 mm connection bolt, right-hand side, - Adjustable dosage dial and pressure gauge - G 3/8 pressure outlet (190 l/min), bottom, with locking screw - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 l/min, UNF 9/16 <p>Connection bolt thread: G 5/8, 90° Bullnose</p>
WM 30857-0000	<p>OXYWAY Fast II, G 5/8, 90° Bullnose Oxygen Pressure Reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G 5/8, 90° Bullnose

	<ul style="list-style-type: none"> - 30 mm connection bolt, right-hand side, - Adjustable dosage dial and pressure gauge - G 3/8 pressure outlet (190 l/min), bottom, with locking screw - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 l/min, UNF 9/16
WM 31891-0000	<p>OXYWAY Fast II, G 3/4 Oxygen pressure reducer,</p> <ul style="list-style-type: none"> - High-pressure manual connection G 3/4 - 30 mm connection bolt, left-hand side - Dosage dial and pressure gauge - G 3/8 pressure outlet G 3/8, bottom (approx. 190 L/min), with locking screw - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16, up 180° turned
WM 31881-0000	<p>OXYWAY Fast II, G 3/4 with flat gauge Oxygen pressure reducer,</p> <ul style="list-style-type: none"> - High-pressure manual connection G 3/4 - 30 mm connection bolt, right-hand side - Dosage dial and pressure gauge (flat, zero-position 3 o'clock) - G 3/8 pressure outlet G 3/8, bottom (approx. 190 L/min), with locking screw - Flow outlet, dosages can be set to 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16
WM 31871-0000	<p>OXYWAY Fast II, G 3/4 with flat gauge Oxygen pressure reducer,</p> <ul style="list-style-type: none"> - High-pressure manual connection G 3/4 - 50 mm connection bolt, right-hand side - Dosage dial and pressure gauge (flat, zero-position 3 o'clock) - G 3/8 pressure outlet G 3/8, bottom - Flow outlet, dosages can be set to 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16, turned 90°, left
WM 33931-0000	<p>OXYWAY Fast II, G 3/4 Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G 3/4 - 30 mm connection bolt, right-hand side - Dosage dial and flat pressure gauge - Pressure outlet: DIN coupling (150 l/min), bottom - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 l/min, UNF 9/16, bottom
WM 33951-0000	<p>OXYWAY Fast II, G 3/4 with flat manometer Oxygen pressure reducer</p> <ul style="list-style-type: none"> - High-pressure manual connection G 3/4 - 50 mm connection bolt, right-hand side - Dosage dial and flat pressure gauge - Pressure outlet: DIN coupling (150 l/min), bottom - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 l/min, UNF 9/16
WM 33791-0000	<p>OXYWAY Fast II, G 3/4 with flat gauge Oxygen pressure reducer,</p> <ul style="list-style-type: none"> - High-pressure manual connection G 3/4 - 30 mm connection bolt, left-hand side - Dosage dial and pressure gauge (flat, zero-position 12 o'clock) - G 3/8 pressure outlet G 3/8, bottom, free flow (approx. 190 L/min), with locking screw - Flow outlet, dosages can be set to 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16, turned 180°, up
WM 31901	<p>OXYWAY Fast II, G 3/4 Oxygen pressure reducer</p>

	<ul style="list-style-type: none"> - High-pressure manual connection G 3/4 - 30 mm connection bolt, left-hand side - Dosage dial and pressure gauge - Pressure outlet: DIN coupling (150 l/min), bottom - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 l/min, UNF 9/16
WM 33981-0000	<p>OXYWAY Fast II, G 3/4 Oxygen pressure reducer,</p> <ul style="list-style-type: none"> - High-pressure manual connection G 3/4 - 20 mm connection bolt, left-hand side - Dosage dial and pressure gauge - G 3/8 pressure outlet G 3/8, bottom, free flow (approx. 190 L/min), with locking screw - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16, up 180° turned
WM 31811-0000	DM_R
WM 33971-0000	<p>OXYWAY Fast II, G 3/4 with flat gauge Oxygen pressure reducer,</p> <ul style="list-style-type: none"> - High-pressure manual connection G 3/4 - 30 mm connection bolt, right-hand side - Dosage dial and pressure gauge (flat, zero-position 3 o'clock) - G 3/8 pressure outlet G 3/8, bottom, free flow (approx. 190 L/min), with locking screw - Flow outlet, dosages can be set to 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16, turned 90°, left
WM 31892-0000	<p>OXYWAY Fast II, G 5/8, 60° Bullnose Oxygen pressure reducer,</p> <ul style="list-style-type: none"> - High-pressure manual connection G 5/8, 60° Bullnose - 30 mm connection bolt, left-hand side - Dosage dial and pressure gauge - G 3/8 pressure outlet G 3/8, bottom (approx. 190 L/min), with locking screw - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16 turned 180°, up
WM 31893-0000	<p>OXYWAY Fast II, G 5/8, 72° Bullnose Oxygen pressure reducer,</p> <ul style="list-style-type: none"> - High-pressure manual connection G 5/8, 72° Bullnose - 30 mm connection bolt, left-hand side - Dosage dial and pressure gauge - G 3/8 pressure outlet G 3/8, bottom (approx. 190 L/min), with locking screw - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16 up 180° turned
WM 31894-0000	<p>"OXYWAY Fast II, W21,8 x 1/14" Oxygen pressure reducer,</p> <ul style="list-style-type: none"> - High-pressure manual connection W21,8 x 1/14" - 30 mm connection bolt, left-hand side - Dosage dial and pressure gauge - G 3/8 pressure outlet G 3/8, bottom (approx. 190 L/min), with locking screw - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16", up 180° turned
WM 31895-0000	<p>OXYWAY Fast II, CGA 540, NGO Oxygen pressure reducer,</p> <ul style="list-style-type: none"> - High-pressure manual connection CGA 540, NGO - 30 mm connection bolt, left-hand side - Dosage dial and pressure gauge

- G 3/8 pressure outlet G 3/8, bottom (approx. 190 L/min), with locking screw
- Flow outlet, dosages can be set to 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16", up 180° turned

- WM 31896-0000 OXYWAY Fast II, PIN-Index
Oxygen pressure reducer,
- High-pressure PIN-Index manual connection
 - 30 mm connection bolt, left-hand side
 - Dosage dial and pressure gauge
 - G 3/8 pressure outlet G 3/8, bottom (approx. 190 L/min), with locking screw
 - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16", up 180° turned

- WM 31897-0000 OXYWAY Fast II, G 5/8, 90° Bullnose
Oxygen pressure reducer,
- High-pressure manual connection G 5/8, 90° Bullnose
 - 30 mm connection bolt, left-hand side
 - Dosage dial and pressure gauge
 - G 3/8 pressure outlet G 3/8, bottom (approx. 190 L/min), with locking screw
 - Flow outlet, adjustable in levels of 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min, UNF 9/16", up 180° turned

OXYWAY Fast III

Connection bolt thread: G 3/4

- WM 30601-0000 OXYWAY Fast III, G 3/4
Oxygen pressure reducer
- High-pressure manual connection G 3/4
 - 30 mm connection bolt, right-hand side
 - Dosage dial and pressure gauge
 - Flow outlet, bottom, UNF 9/16
 - Dosages in levels 0,5, 1, 1,5, 2, 2,5, 3, 4, 5, 6 l/min

- WM 30631-0000 OXYWAY Fast III, G 3/4
Oxygen pressure reducer
- High-pressure manual connection G 3/4
 - 80 mm connection bolt, right-hand side
 - Dosage dial and pressure gauge
 - Flow outlet, bottom, UNF 9/16
 - Dosages in levels 0,5, 1, 1,5, 2, 2,5, 3, 4, 5, 6 l/min

- WM 30611-0000 OXYWAY Fast III, G 3/4
Oxygen pressure reducer
- High-pressure manual connection G 3/4
 - 30 mm connection bolt, right-hand side
 - Dosage dial and pressure gauge, flat
 - Flow outlet, bottom, UNF 9/16
 - Dosages in levels 0,5, 1, 1,5, 2, 2,5, 3, 4, 5, 6 l/min

Connection bolt thread: G 5/8,60° Bullnose

- WM 30602-0000 OXYWAY Fast III, G 5/8, 60° Bullnose
Oxygen pressure reducer
- High-pressure manual connection G 5/8, 60° Bullnose
 - 30 mm connection bolt, right-hand side
 - Dosage dial and pressure gauge
 - Flow outlet, bottom, UNF 9/16
 - Dosages in levels 0,5, 1, 1,5, 2, 2,5, 3, 4, 5, 6 l/min

Connection bolt thread: G 5/8,

72° Bullnose

- WM 30623-0000 OXYWAY Fast III, G 5/8, 72° Bullnose
Oxygen pressure reducer
- High-pressure manual connection G 5/8, 72° Bullnose
 - 50 mm connection bolt, right-hand side
 - Dosage dial and pressure gauge
 - Flow outlet, bottom, UNF 9/16
 - Dosages in levels 0,5, 1, 1,5, 2, 2,5, 3, 4, 5, 6 l/min

Connection bolt thread: W21,8 x 1/14"

- WM 30604-0000 "OXYWAY Fast III, W21,8 x 1/14"
Oxygen pressure reducer
- High-pressure manual connection W21,8 x 1/14"
 - 30 mm connection bolt, right-hand side
 - Dosage dial and pressure gauge
 - Flow outlet, bottom, UNF 9/16
 - Dosages in levels 0,5, 1, 1,5, 2, 2,5, 3, 4, 5, 6 l/min

OXYWAY Flow Dosage

OXYWAY Click

- WM 31030-0000 OXYWAY Click
with DIN connector (DIN 13260),
Flow dosage setting knob for central gas system
in emergency vehicles, clicks into place at
settings for 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min,
incl. one (1) connection nozzle with swivel nut,
straight, Ø 6 mm, UNF 9/16 (WM 31122)

- WM 31040-0000 OXYWAY Click
with AGA connector,
Flow dosage setting knob for central gas facility
in emergency vehicles, clicks into place at
settings for 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min,
incl. one (1) connection nozzle with swivel nut,
straight, Ø 6 mm, UNF 9/16 (WM 31122)

- WM 31050-0000 OXYWAY Click
with SIS connector (AS 2896),
Flow dosage setting knob for central gas system,
clicks into place at
settings for 1, 2, 3, 4, 5, 6, 9, 12, 15 L/min,
incl. one (1) connection nozzle with swivel nut,
straight, Ø 6 mm, UNF 9/16 (WM 31122)

Accessories and Replacement Parts for all pressure reducers

- WM 31122 Connection nozzle with swivel nut, straight,
Ø 6 mm, UNF 9/16
- WM 31130 Angled connection nozzle with swivel nut,
90°, Ø 6 mm, UNF 9/16
- WM 31140 Angled connection nozzle with wing nut,
90°, Ø 6 mm, UNF 9/16
- WM 1415 Connection nozzle with swivel nut, Ø 6,7 mm,
G 3/8 (for WM 1411)
- WM 1470 Straight connection nozzle with swivel nut,
Ø 5.2 mm, G 3/8
- WM 1475 Angled connection nozzle with wing nut,
Ø 5.2 mm, G 3/8

WM 22319	Angled connection nozzle with swivel nut, angled 90°, Ø 5.2 mm, G 3/8
WM 1231	Blind cap for G 3/8 male threads
WM 29148	Blind cap for UNF 9/16 male threads
WM 1463	Locking screw with G 3/8 male thread
WM 30905	Filter screw Change every 2 years
WM 30008	Torx screw for housing of pressure reducer gauge with base (to 06/2001)
WM 30009	Torx screw for housing of pressure reducer gauge with base (from 06/2001)
WM 1145	O-Ring, 12.3 x 2.4
WM 1148	Set O-rings for connection G 3/4
WM 31111	Gauge cap to protect from impact and falls for gauge with wrench size 14 mm
WM 31024	Gauge cap to protect from impact and falls for flat contents gauge 0 - 315 bar
WM 12825	Bubble humidifier with UNF 9/16" thread

Consumables for inhalation

WM 1429	Inhalation mask, plastic, with hose The article contains the following substances included on the candidate list according to article 59 (1, 10) of Regulation (EC) No 1907/2006 (REACH) in a concentration above 0.1 % weight by DEHP (CAS No. 117-81-7, EC No. 204-211-0) The packaging materials contain no substances included on the candidate list in a concentration above 0.1 % weight by
WM 1925	Nasal cannula for adults, two-pronged with harness and 2,1 m connection hose (nasal cannula) The article contains the following substances included on the candidate list according to article 59 (1, 10) of Regulation (EC) No 1907/2006 (REACH) in a concentration above 0.1 % weight by weight: DEHP (CAS No. 117-81-7, EC No. 204-211-0) The packaging materials contain no substances included on the candidate list in a concentration above 0.1 % weight by weight.
WM 1929	Nasal cannula for adults, two-pronged with harness and 1,5 m connection hose (nasal cannula) The article contains the following substances included on the candidate list according to article 59 (1, 10) of Regulation (EC) No 1907/2006 (REACH) in a concentration above 0.1 % weight by weight: DEHP (CAS No. 117-81-7, EC No. 204-211-0) The packaging materials contain no substances included on the candidate list in a concentration above 0.1 % weight by weight.

Oxygen Control Panels

Oxygen administrators

WM 3714	Multiple oxygen administrator with two (2) outlets,
---------	---

	reverse screw-on type Inlet: one (1) inlet G 3/8 with male thread, top Outlet: two (2) outlets G 3/8 with male thread, free
WM 3730	Multiple oxygen administrator with 3 outlets Inlet: one (1) inlet wing nut G 3/8, top Outlet: three (3) outlets G 3/8 with male thread, free
WM 3834	Multiple oxygen administrator with 3 outlets Inlet: one (1) G3/8 Outlet: three (3) UNF 9/16, 10 L/min including three (3) angled connection nozzle
WM 3740	Multiple oxygen administrator with four (4) outlets Inlet: one (1) G 3/8 inlet wing nut, top Outlet: four (4) outlets G 3/8 with male thread, free
WM 3880-0000	OXYWAY Click multiple oxygen administrator with four (4) dosage setting stops for outlets Outlet: Four (4) UNF 9/16 male thread, 0-15 l/min, dosage setting stops, incl. four (4) connection nozzle UNF 9/16 (WM 31122) Inlet / outlet: Three (3) G 3/8 male thread, bottom, top, right
WM 3883	OXYWAY Click multiple oxygen administrator with four (4) dosage setting stops for outlets Outlet: Four (4) UNF 9/16 male thread, 0-15 l/min, dosage setting stops, incl. four (4) angled connection nozzle UNF 9/16 (WM 31130) Inlet / outlet: One (1) G 3/8 male thread, top Two (2) adapter on coupling, t
WM 3899	OXYWAY Click multiple oxygen administrator with four (4) dosage setting stops for outlets Outlet: Four (4) UNF 9/16 male thread, 0-15 l/min, dosage setting stops, incl. four (4) angled connection nozzle UNF 9/16 (WM 31130) Inlet / outlet: One (1) G 3/8 male thread, One (1) T-shaped distributor
WM 3804-0000	OXYWAY Click multiple oxygen administrator with four (4) dosage setting stops for outlets and Multiport, with two adapters for standard device rails Outlet: Four (4) UNF 9/16 male thread, 0-15 l/min, dosage setting stops, incl. four (4) angled connection nozzles UNF 9/16 (WM 31130) Inlet / outlet: 1 x G 3/8 male thread

Accessories and Replacement Parts for Oxygen Control Panels

O2 plugs and adapters

WM 28627	Connector (BS 5682) for central gas system outlet, straight (accessories required for assembly: WM 2078)
----------	--

WM 22345	Connector (AGA) for central gas system outlet (required for assembly: ear clip WM 2078)
WM 22252	Angled bayonet connector (AFNOR) for central gas system outlet
WM 22296	Straight O2 plug type bayonet (AFNOR) for central gas facility outlet
WM 22929	SIS screw connector for central gas facilities in compliance with Australian standard (AS 2896), with hose connection nozzle
WM 2063	O2 connector (UNI 9507) for central gas system outlet
WM 4847	Straight O2 connector (EMS), Ø 5,5, chrome-plated
WM 29141	Quick release coupler type Rectus, series 21, with G 1/8 external thread (e.g. for MEDUMAT Transport, Modules)
WM 29149	Quick release coupler type Rectus, series 21, with G 3/8 internal thread (e.g. for OXYWAY pressure reducers)
WM 29150	Quick release coupler type Rectus, series 21, with connection nozzle Ø 6 mm
WM 29151	O2 connector type Rectus, series 21, with connection nozzle Ø 6 mm
WM 29131	O1 Connector type OHMEDA, green, with straight connection nozzle (required for assembly: ear clip WM 2078)
WM 29126	Adapter for quick release coupler type Ohmeda, green, with DISS screw connector (UNF 9/16 male thread)
WM 4361	Walther coupling with hose connection
WM 2057	O2 connector (DIN 13260) for central gas system outlet (accessories required for assembly: WM 2078)
WM 3799	Plug connector for quick-release coupler type Walther, angled 90°: - Seal nipple type Walther, nominal diameter 4 mm, - hose connection Ø 6 mm (required accessories for assembly: WM 2078)
WM 2062	O2 connector (Carbamed) for central gas system outlet
WM 4364	Seal nipple for oxygen cap (type Walther), straight (required accessories for assembly: WM 2078)
WM 3881	Adapter with G 3/8 internal screw thread and DIN quick release coupler
WM 3949	Adapter for plug (DIN 13260) with aluminium handle
WM 3950-0000	DIN adapter from swivel nut G 3/8 female thread with straight connection nozzle to O2 quick- release coupler type DIN 13260, with 20 cm hose connection
WM 1485	Adapter with G 3/8 internal screw thread and quick release coupler type Walther
WM 15554	Set adapter G 3/8 internal screw thread on NIST-external screw thread for MEDUMATs, modules and OXYWAY pressure reducers
WM 22584	Adapter G 3/8 internal screw thread on NIST- external screw thread for MEDUMATs, modules

	and OXYWAY pressure reducers
WM 22906	Adapter, G3/8 - G3/8, chrome-plated
WM 15575	Adapter set G 3/8 to M 12 to connect Weinmann pressure reducer to Dräger pressure hose
WM 13755	"Adapter for G3/8" internal thread to UNF 9/16" external thread "
WM 22217	Adapter M12X1 - G3/8, chrome-plated
WM 1412	Screw socket with male thread G 3/8
WM 1145/31	Sealing washer, 7 x 14 x 2, for 3/8 threaded connection
WM 3797	Holder for multiple oxygen administrator
WM 22905	Adapter on coupling, type bayonet (AFNOR)
WM 22910	Angular adapter on coupling, type bayonet (AFNOR)
WM 30048	Angular adapter with DIN 13260 for G3/8" thread and DIN 13260 coupling
WM 1490	Adapter with G 3/8 internal screw thread and quick release coupler for British Standard 5682 (Schrader)

Connection Nozzles

WM 31122	Connection nozzle with swivel nut, straight, Ø 6 mm, UNF 9/16
WM 31130	Angled connection nozzle with swivel nut, 90°, Ø 6 mm, UNF 9/16
WM 29140	Angled connection nozzle with swivel nut, pulling protection, 90°, Ø 6 mm, UNF 9/16
WM 31140	Angled connection nozzle with wing nut, 90°, Ø 6 mm, UNF 9/16
WM 1415	Connection nozzle with swivel nut, Ø 6,7 mm, G 3/8 (for WM 1411)
WM 1470	Straight connection nozzle with swivel nut, Ø 5.2 mm, G 3/8
WM 1475	Angled connection nozzle with wing nut, Ø 5.2 mm, G 3/8
WM 22319	Angled connection nozzle with swivel nut, angled 90°, Ø 5.2 mm, G 3/8
WM 22691	Screw socket NIST

Distributors

WM 3816	Shut-off valve with one outlet Inlet: one (1) inlet wing nut G3/8, top Outlet: one (1) outlet G3/8 with male thread, secretion, 13 L/min
WM 3911	T-shaped multiple oxygen administrator, two outlets G 3/8, one inlet wing nut
WM 1289	T-shaped distributor two (2) outlets G 3/8, one (1) inlet wing nut
WM 3910	T-shaped distributor, one (1) outlet G 3/8, one (1) inlet G3/8,

- one (1) inlet cap (type Walther)
- WM 8220 T-shaped distributor with check valve
one (1) outlet G 3/8,
one (1) inlet cap (type Walther)
- WM 1472 Connection nozzle with wing nut G 3/8 and two outlets Ø 5.2 mm

Pressure hoses

Alternative hosing upon request.

- WM 2081 "Low pressure hose, WHITE, medical PVC
(DEHP free), wall thickness 3 mm, inner diameter
5 mm, marking: ""Med. PVC, 1000kPa, 10bar,
WEINMANN""
"
- WM 2069 "Low pressure hose, BLACK, medical PVC
(DEHP free), wall thickness 3 mm, inner diameter
5 mm, marking: ""Med. PVC, 1000kPa, 10bar,
WEINMANN""
"
- WM 2068 "Low pressure hose, GREEN, medical PVC
(DEHP free), wall thickness 3 mm, inner diameter
5 mm, marking: ""Med. PVC, 1000kPa, 10bar,
WEINMANN""
"
- WM 36300-0000 Pressure hose with
2x swivel nut G 3/8 female thread with angled
connection nozzle
Length: 1m, white
- WM 36301-0000 "Pressure hose with
1x swivel nut G 3/8"" female thread with straight
connection nozzle,
1x swivel nut G 3/8"" female thread with angled
connection nozzle
"
Length: 1m, white
- WM 36302-0000 Pressure hose with
1x swivel nut G 3/8 inner thread with angled
connection nozzle,
1x angled O2 connector for quick release coupler
type Walther
Length: 1m, white
- WM 36303-0000 "Pressure hose with
1x swivel nut G 3/8"" inner thread with angled
connection nozzle,
1x straight O2 connector for quick release coupler
type Walther
"
Length: 1m, white
- WM 36304-0000 "Pressure hose with
1x swivel nut G 3/8"" inner thread with angled
connection nozzle,
1x Walther coupling
"
Length: 1m, white
- WM 36305-0000 "Pressure hose with
1x swivel nut G 3/8"" inner thread with angled
connection nozzle,
1x angled O2 connector for quick release coupler
standard DIN 13260
"
Length: 1m, white

WM 36306-0000	"Pressure hose with 1x swivel nut G 3/8"" inner thread with angled connection nozzle, 1x angled O2 connector for quick release coupler standard NF-S 90-116 (AFNOR, Bayonett) " Length: 1m, white
WM 36307-0000	"Pressure hose with 1x swivel nut G 3/8"" inner thread with angled connection nozzle, 1x straight O2 connector for quick release coupler standard NF-S 90-116 (AFNOR, Bayonett) " Length: 1m, white
WM 36308-0000	"Pressure hose with 1x swivel nut G 3/8"" inner thread with angled connection nozzle, 1x straight O2 connector for quick release coupler standard SS 87 524 30 (AGA) " Length: 1m, white
WM 36309-0000	"Pressure hose with 1x swivel nut G 3/8"" inner thread with angled connection nozzle, 1x straight O2 connector for quick release coupler standard AS 2896 " Length: 1m, white
WM 36310-0000	"Pressure hose with 1x swivel nut G 3/8"" female thread with angled connection nozzle, 1x swivel nut UNF 9/16"" female thread with angled connection nozzle " Length: 1m, white
WM 36311-0000	"Pressure hose with 1x swivel nut G 3/8"" inner thread with angled connection nozzle, 1x angled O2 connector for quick release coupler type Carbamed " Length: 1m, white
WM 36312-0000	"Pressure hose with 1x swivel nut G 3/8"" inner thread with angled connection nozzle, 1x straight O2 connector for quick release coupler standard BS 5682 (Schrader, BOC) " Length: 1m, white
WM 36313-0000	"Pressure hose with 1x swivel nut G 3/8"" inner thread with angled connection nozzle, 1x straight O2 connector for quick release coupler type Ohmeda " Length: 1m, white
WM 36314-0000	"Pressure hose with 1x swivel nut G 3/8"" inner thread with angled connection nozzle, 1x straight O2 connector for quick release coupler

	standard UNI 9507 "
	Length: 1m, white
WM 36315-0000	"Pressure hose with 2x swivel nut G 3/8"" female thread with straight connection nozzle "
	Length: 1m, white
WM 36316-0000	"Pressure hose with 1x swivel nut G 3/8"" inner thread with straight connection nozzle, 1x angled O2 connector for quick release coupler type Walther "
	Length: 1m, white
WM 36317-0000	Pressure hose with 1x Walther coupling, 1x angled O2 connector for quick release coupler type Walther Length: 1m, white
WM 36318-0000	"Pressure hose with 1x swivel nut G 3/8"" inner thread with straight connection nozzle, 1x angled O2 connector for quick release coupler standard DIN 13260 "
	Length: 1m, white
WM 36319-0000	Pressure hose with 1x angled O2 connector for quick release coupler standard DIN 13260 1x angled O2 connector for quick release coupler type Walther Length: 1m, white
WM 36320-0000	"Pressure hose with 1x swivel nut G 3/8"" inner thread with straight connection nozzle, 1x angled O2 connector for quick release coupler standard NF-S 90-116 (AFNOR, Bayonett) "
	Length: 1m, white
WM 36321-0000	Pressure hose with 1x angled O2 connector for quick release coupler type Walther, 1x angled O2 connector for quick release coupler standard NF-S 90-116 (AFNOR, Bayonett) Length: 1m, white
WM 36322-0000	"Pressure hose with 1x swivel nut G 3/8"" inner thread with straight connection nozzle, 1x angled O2 connector for quick release coupler standard SS 87 524 30 (AGA) "
	Length: 1m, white
WM 36323-0000	Pressure hose with 1x angled O2 connector for quick release coupler type Walther, 1x angled O2 connector for quick release coupler standard SS 87 524 30 (AGA) Length: 1m, white
WM 36324-0000	"Pressure hose with

	1x angled O2 connector for quick release coupler type Walther, 1x swivel nut UNF 9/16" female thread with angled connection nozzle "
	Length: 1m, white
WM 36325-0000	"Pressure hose with 1x straight O2 connector for quick release coupler type Walther, 1x swivel nut UNF 9/16" female thread with angled connection nozzle "
	Length: 1m, white
WM 36326-0000	Pressure hose with 1x straight O2 connector for quick release coupler standard NF-S 90-116 (AFNOR, Bayonett), 1x angled O2 connector for quick release coupler standard NF-S 90-116 (AFNOR, Bayonett) Length: 1m, white
WM 36327-0000	Pressure hose with 1x straight O2 connector for quick release coupler type Ohmeda, 1x straight O2 connector for quick release coupler type Walther Length: 1m, white
WM 36328-0000	"Pressure hose with 1x swivel nut G 3/8" inner thread with angled connection nozzle, 1x straight O2 connector for quick release coupler type Rectus 21 (Medprom) "
	Length: 1m, white
WM 36329	"Pressure hose with 1x swivel nut G 3/8" inner thread with angled connection nozzle, 1x angled O2 connector for quick release coupler standard DIN 13260 1x quick coupling Type Walther "
WM 36331-0000	Pressure hose with 1x angled O2 connector for quick release coupler type Walther 1x straight O2 connector for quick release coupler type BS 5682

Oxygen Cylinders

Further connection threads are available upon request

Oxygen Cylinders 0.3 L

WM 1812	0.3L oxygen cylinder, filled, G 3/4, max. pressure 200 bar
WM 1812-02	0.3L oxygen cylinder, empty, G 3/4, max. pressure 200 bar
WM 1812/2-02	0.3L oxygen cylinder, empty, G 5/8 Bullnose 90°, white, max. pressure 200 bar
WM 1812/12-02	0.3L oxygen cylinder, empty, G 5/8 Bullnose 72°, white, max. pressure 200 bar
WM 1812/5-02	0.3L oxygen cylinder, empty, G 5/8 Bullnose 60°, black, O2med, max. pressure 200 bar

WM 1812/06-02 0,3L oxygen cylinder, empty, NGO, green, O2med, with safety valve 250 bar, max. pressure 200 bar

Oxygen Cylinders 0.8 L

WM 1818 0.8L oxygen cylinder, filled, G 3/4, max. pressure 200 bar

WM 1818-02 0.8L oxygen cylinder, empty, G 3/4, max. pressure 200 bar

WM 1818/005-02 0.8L oxygen cylinder, empty, PIN-Index, black, O2med, max. pressure 200 bar

WM 1818/006-02 0.8 L Oxygen cylinder, empty, 200 bar, PIN-Index, green, O2med

WM 1818/008-02 0.8L oxygen cylinder, empty, PIN-Index, white, O2med, max. pressure 200 bar

WM 1818/2-02 0.8L oxygen cylinder, empty, G 5/8 Bullnose 90°, white, max. pressure 200 bar

WM 1818/4-02 0.8L oxygen cylinder, empty, W21,8x1/14, white, O2med, max. pressure 200 bar

WM 1818/5-02 0.8L oxygen cylinder, empty, G 5/8 Bullnose 60°, black, O2med, max. pressure 200 bar

WM 1818/6-02 0.8L oxygen cylinder, empty, NGO, green, O2med, max. pressure 200 bar

WM 1818/06-02 0,8L oxygen cylinder, empty, NGO, green, O2med, with safety valve 250 bar, max. pressure 200 bar

WM 1818/7-02 0.8L oxygen cylinder, empty, G 5/8 Bullnose 60°, white, max. pressure 200 bar

WM 1818/8-02 0.8L oxygen cylinder, empty, G 5/8 Bullnose 60°, white, O2med, max. pressure 200 bar

WM 1818/10-02 0.8L oxygen cylinder, empty, W21,8x1/14, white, max. pressure 200 bar

WM 1818/11-02 0.8L oxygen cylinder, empty, W21,8x1/14, white, with collar, max. pressure 200 bar

WM 1818/12-02 0.8L oxygen cylinder, empty, G 5/8 Bullnose 72°, white, max. pressure 200 bar

WM 1818/14-02 0.8L oxygen cylinder, empty, G 5/8 Bullnose 72°, white, O2med, max. pressure 200 bar

WM 1818/19-02 0.8L oxygen cylinder, empty, G 3/4, blue, O2med, max. pressure 200 bar

Oxygen Cylinders 2 L

WM 1822 2L oxygen cylinder, filled, G 3/4, max. pressure 200 bar

WM 1822-02 2L oxygen cylinder, empty, G 3/4, max. pressure 200 bar

WM 1822/005-02 2L oxygen cylinder, empty, PIN-Index, black, O2med, max. pressure 200 bar

WM 1822/006-02 2L Oxygen cylinder, empty, 200 bar, PIN-Index, green, O2med

WM 1822/008-02 2L oxygen cylinder, empty, PIN-Index, white, O2med, max. pressure 200 bar

WM 1822/2-02 2L oxygen cylinder, empty, G 5/8 Bullnose 90°, white, max. pressure 200 bar

WM 1822/4-02 2L oxygen cylinder, empty, W21,8x1/14, white, O2med, max. pressure 200 bar

WM 1822/5-02	2L oxygen cylinder, empty, G 5/8 Bullnose 60°, black, O2med, max. pressure 200 bar
WM 1822/6-02	2L oxygen cylinder, empty, NGO, green, O2med, max. pressure 200 bar
WM 1822/06-02	2L oxygen cylinder, empty, NGO, green, O2med, with safety valve 250 bar, max. pressure 200 bar
WM 1822/7-02	2L oxygen cylinder, empty, G 5/8 Bullnose 60°, white, max. pressure 200 bar
WM 1822/8-02	2L oxygen cylinder, empty, G 5/8 Bullnose 60°, white, O2med, max. pressure 200 bar
WM 1822/10-02	2L oxygen cylinder, empty, W21,8x1/14, white, max. pressure 200 bar
WM 1822/11-02	2L oxygen cylinder, empty, W21,8x1/14, white, with collar, max. pressure 200 bar
WM 1822/12-02	2L oxygen cylinder, empty, G 5/8 Bullnose 72°, white, max. pressure 200 bar
WM 1822/14-02	2L oxygen cylinder, empty, G 5/8 Bullnose 72°, white, O2med, max. pressure 200 bar
WM 1822/18-02	2L oxygen cylinder, empty, W21,8x1/14, green, O2med, max. pressure 200 bar
WM 1822/19-02	2L oxygen cylinder, empty, G 3/4, blue, O2med, max. pressure 200 bar
WM 1822/20-02	2L oxygen cylinder, empty, NF E 29-650/F, white, max. pressure 200 bar

Oxygen Cylinders 2 L Aluminum

WM 1814	2 l oxygen lightweight-cylinder, filled, G 3/4, aluminium, max. pressure 200 bar
WM 1814-02	2 l oxygen lightweight-cylinder, empty, G 3/4, aluminium, max. pressure 200 bar
WM 1814/6-02	2L oxygen lightweight-cylinder, empty, NGO, aluminium, green, O2-Med, max. pressure 200 bar
WM 1814/06-02	2L oxygen lightweight-cylinder, empty, NGO, aluminium, green, O2med, with safety valve 250 bar, max. pressure 200 bar
WM 1814/18-02	2L oxygen lightweight-cylinder, empty, W21,8x1/14, green, aluminum, O2med, max. pressure 200 bar

Oxygen Cylinders 5 L

WM 1805	5L oxygen cylinder, filled, G 3/4, max. pressure 200 bar
WM 1805-02	5L oxygen cylinder, empty, G 3/4, max. pressure 200 bar
WM 1805/6-02	5L oxygen cylinder, empty, NGO, green, O2med, max. pressure 200 bar
WM 1805/4-02	5L oxygen cylinder, empty, W21,8x1/14, white, O2 med, max. pressure 200 bar
WM 1805/8-02	5L oxygen cylinder, empty, G 5/8 Bullnose 60°, white, O2med, max. pressure 200 bar
WM 1805/10-02	5L oxygen cylinder, empty, W21,8x1/14, white, max. pressure 200 bar
WM 1805/11-02	5L oxygen cylinder, empty, W21,8x1/14, white, with collar, max. pressure 200 bar

Oxygen Cylinders 10 L

WM 1825	10L oxygen cylinder, filled, G 3/4, with safety cap and base ring, max. pressure 200 bar
WM 1825-02	10L oxygen cylinder, empty, G 3/4, with safety cap and base ring, max. pressure 200 bar
WM 1825/4-02	10L oxygen cylinder, empty, W21,8x1/14, white, with safety cap and integrated base, max. pressure 200 bar
WM 1825/5-02	10L oxygen cylinder, empty, G 5/8 Bullnose 60°, black, O2med, with safety cap and integrated base, max. pressure 200 bar
WM 1825/6-02	10L oxygen cylinder, empty, NGO, green, O2med, with safety cap and integrated base, max. pressure 200 bar
WM 1825/7-02	10L oxygen cylinder, empty, G 5/8 Bullnose 60°, white, with safety cap and integrated base, max. pressure 200 bar
WM 1825/8-02	10L oxygen cylinder, empty, G 5/8 Bullnose 60°, white, O2med, filled, with safety cap and integrated base, max. pressure 200 bar
WM 1825/10-02	10L oxygen cylinder, empty, W21,8x1/14, white, with safety cap and integrated base, max. pressure 200 bar
WM 1825/11-02	10L oxygen cylinder, empty, W21,8x1/14, white, with safety cap and base ring, max. pressure 200 bar
WM 1825/12-02	10L oxygen cylinder, empty, G 5/8 Bullnose 72°, white, with safety cap and integrated base, max. pressure 200 bar
WM 1825/14-02	10L oxygen cylinder, empty, G 5/8 Bullnose 72°, white, O2med, with safety cap and integrated base, max. pressure 200 bar
WM 1825/19-02	10L oxygen cylinder, empty, G 3/4, blue, O2med, with safety cap and integrated base, max. pressure 200 bar
WM 1825/20-02	10L oxygen cylinder, empty, NF E 29-650/F, white, with safety cap and integrated base, max. pressure 200 bar

Accessories for Oxygen Cylinders

Adapters for Oxygen Cylinders

WM 2718	Adapter G3/4, chrome-plated (foreign 5/8, 72° internal thread) for: Algeria, Belgium, France, Holland, Portugal, Spain
WM 2732	Adapter G3/4 (foreign W 21,8) for: Denmark, Finland, Italy, Ex-Yugoslavia, Luxembourg, Norway, Syria, Czech Republic, Turkey (Please order WM 2831 separately)
WM 17815	Adapter Set G3/4 (foreign W 21,8) for: Denmark, Finland, Italy, Ex-Yugoslavia, Luxembourg, Norway, Syria, Czech Republic, Turkey
WM 2735	Adapter G3/4 (foreign 5/8 BSP) for: Arabia, Egypt, England, Greece, India, Indonesia,

- Pakistan, South Africa
- WM 2714 Adapter G3/4 (foreign 5/8 BSP) for:
Arabia, Egypt, England, Greece, India, Indonesia,
Pakistan, South Africa
- WM 2751 Adapter G3/4 (foreign 5/8 external thread) for:
Algeria, Belgium, France, Holland, Portugal,
Spain
- WM 2753 Adapter G3/4 (foreign W21,8) for:
Denmark, Finland, Italy, Ex-Yugoslavia,
Luxembourg, Norway, Syria, Czech Republic,
Turkey
- WM 2769 Adapter G3/4 (foreign NGO) for:
Ecuador, Japan, Canada, Korea, Peru,
San Salvador, Thailand, USA

Cylinder Brackets

- WM 1885 Set of 2 cylinder brackets with console,
rubber profile, set of two for stationary O2 facility
in buildings or vehicles for 2.0 and
3.0-L O2 cylinders,
100 mm Ø (complies with EN 1789 norm)
- WM 1889 Cylinder bracket with console, rubber profile,
set of two for stationary O2 facility in buildings
or vehicles for 5.0, 10.0 or 11.0-L O2 cylinder,
140 mm Ø (complies with EN 1789 norm)
- WM 15071 Set of straining screws for cylinder brackets
Ø 83, Ø 100, Ø 115, with console and rapid
action tension jack
- WM 15072 Set of straining screws for cylinder brackets
Ø 140, with console and rapid action tension jack
- WM 15073 Set of straining screws for cylinder brackets
Ø 100, Ø 115, without console and Ø 115 with
console
- WM 15074 Set of straining screws for cylinder brackets
Ø 83, without console

OXYBAG - bags for oxygen bottles

- WM 9802 OXYBAG for 2-Liter oxygen cylinder
(Colour: Blue)
(Max. bottle length incl. pressure reducer: 500 mm)
- WM 9805 OXYBAG for 2 liter to 5 liter oxygen cylinder
(Colour: Blue)
(Max. bottle length incl. pressure reducer: 550 mm)
- WM 9812 OXYBAG for 2 Liter oxygen cylinder
(Colour: Green)
(Max. bottle length incl. pressure reducer: 500 mm)
- WM 9815 OXYBAG for 2 liter to 5 liter oxygen cylinder
(Colour: Green)
(Max. bottle length incl. pressure reducer: 550 mm)

O2 Box Rescue

Please note that the prices for modular systems containing oxygen bottles or pressure reducers may vary if delivered with other connection threads.

- WM 9500-0000 O2 Box Rescue,
Portable, multi-use oxygen system

Box type Zarges, dimensions: 620 x 430 x 250 mm:
Two (2) oxygen cylinders, 5.0-liters x 200 bar, empty
(WM 1805-02)
Two (2) OXYWAY Fix III pressure reducers type Walther (WM 30301)
Four (4) cylinder brackets (WM 1875)
One (1) OXYWAY Click multiple oxygen administrator with four (4) dosage setting stops for outlets (WM 3899)
Four (4) oxygen inhalation masks with connection hose (WM 1429), packed in two (2) accessory bags

Weight including oxygen cylinders: approx. 28 kg
Weight without oxygen cylinders: approx. 6.5 kg

WM 9515-0000

O2 Box Rescue compact
consists of:
One (1) aluminium box, stackable, with hole pattern and foam insert
Outer diameter (WxHxD): 615 x 420 x 270 mm
Inner diameter (BxHxT): 570 x 380 x 220 mm
Weight: approx. 5 kg
Equipped with one carrying handle, cover with two belts and cover holder
Two (2) 5 L oxygen cylinder, empty, G 3/4,
Two (2) OXYWAY Fix III, G 3/4,
One (1) pressure hose, 1 m, white, with plug connector for quick-release coupler type Walther and angled connection nozzle with swivel nut G 3/8,
4x inhalation mask, synthetic, with hose (WM 1429),
2x accessory bag, size M, blue, 250 x 120 x 70 mm (WM 9026),
1x holding plate for device rail adapters (WM 3805),
OXYWAY Click multiple oxygen administrator

WM 9530-0000

O2 Box Rescue Ventilation
consists of:
One (1) aluminium box, stackable, with hole pattern and foam insert
Outer diameter (WxHxD): 600 x 340 x 400 mm
Inner diameter (BxHxT): 550 x 310 x 305 mm
Weight: approx. 5 kg
Equipped with one carrying handle each on the left- and right-hand side, cover with two belts and stainless steel hinges,
One (1) 5 L oxygen cylinder, empty, G 3/4,
One (1) OXYWAY Fix III, G 3/4,
One (1) pressure hose, 38 cm, white, with plug connector for quick-release coupler type Walther and angled connection nozzle with swivel nut G 3/8,
One (1) MEDUMAT Standard a (WM 22800),
One (1) MANUVAC (WM 10800),
Set of disposable suction catheters (WM 15051), consists of three of each of the following:
Ch. 06 / Ø 2 mm,
Ch. 10 / Ø 3 mm,
Ch. 12 / Ø 4 mm .
One (1) secretions nozzle with fingertip control and stopper (WM 2251),

One (1) silicone ventilation mask, size 0
(small child) (WM 11110)
One (1) silicone ventilation mask, size 3 (child)
(WM 11113)
One (1) silicone ventilation mask, size 5 (adult)
(WM 11115)
One (1) inhalation mask with bag and hose
(WM 1422),
Ten (10) fleece cloth, 200 x 300 mm (1x WM
3124),
One (1) set of oropharyngeal tubes (Guedel)
(WM 15483)
consists of one each of the following:
- for adults, Size 3,
- for juveniles, Size 2,
- for children, Size 1
- for small children, Size 0
- for neonates, Size 000
One (1) set of nasalpharyngeal tubes (Wendl)
(WM 15484)
consists of one each of the following:
- ch. 26,
- ch. 30

The article contains the following
substances included on the candidate
list according to article 59 (1, 10) of
Regulation (EC) No 1907/2006 (REACH) in a
concentration above 0.1 % weight by weight:
DEHP (CAS No. 117-81-7, EC No. 204-211-0)
The packaging materials contain no
substances included on the candidate list in a
concentration above 0.1 % weight by weight.

Accessories and Replacement Parts

WM 3899 OXYWAY Click
multiple oxygen administrator with
four (4) dosage setting stops for outlets
Outlet:
Four (4) UNF 9/16 male thread, 0-15 l/min, dosage setting
stops, incl. four (4) angled connection nozzle
UNF 9/16 (WM 31130)
Inlet / outlet:
One (1) G 3/8 male thread,
One (1) T-shaped distributor

WM 3804-0000 OXYWAY Click
multiple oxygen administrator with
four (4) dosage setting stops for outlets
and Multiport,
with two adapters for standard device rails
Outlet:
Four (4) UNF 9/16 male thread, 0-15 l/min,
dosage setting stops,
incl. four (4) angled connection
nozzles UNF 9/16 (WM 31130)
Inlet / outlet:
1 x G 3/8 male thread